




# CS-Pro VET

## Equine Protocol Recommendation

### Protocols

- 1. External Occipital Protuberance**  
Standoff: 20 mm  
EFD: 0.19 mJ/mm<sup>2</sup>
- 2. Cervical Articular Facets**  
Standoff: 30 mm  
EFD: 0.22 mJ/mm<sup>2</sup>
- 3. Medial Meniscus of the Stifle**  
Standoff: 30 mm  
EFD: 0.22 mJ/mm<sup>2</sup>
- 4. Scapulothoracic Joint**  
Standoff: 20 mm  
EFD: 0.19 mJ/mm<sup>2</sup>
- 5. Lateral Collateral Ligament of the Elbow**  
Standoff: 20 mm  
EFD: 0.19 mJ/mm<sup>2</sup>
- 6. Bucked Shins**  
Standoff: 5 mm  
EFD: 0.19 mJ/mm<sup>2</sup>
- 7. Proximal Suspensory Ligament**  
Standoff: 30 mm  
EFD: 0.22 mJ/mm<sup>2</sup>
- 8. Tendon Injuries**  
Standoff: 5 mm  
EFD: 0.19 mJ/mm<sup>2</sup>
- 9. Foot Injuries**  
Standoff: 20 mm  
EFD: 0.19 mJ/mm<sup>2</sup>
- 10. Lower Hock Joint**  
Standoff: 30 mm  
EFD: 0.22 mJ/mm<sup>2</sup>

For more information, contact  
[sales@curative-sound.com](mailto:sales@curative-sound.com)



### General Recommendations

Treat each area with 3000 pulses at 10 hertz.

Conduct treatment every 3-10 days for a total of 3-6 treatments.

These recommendations are only general guidelines and can be deviated from at veterinary discretion.

### Other Treatments

#### Superficial Wounds

Standoff: 5 mm  
EFD: 0.19 mJ/mm<sup>2</sup>

#### Splints

Standoff: 20 mm  
EFD: 0.19 mJ/mm<sup>2</sup>