




# CS-Pro VET

30mm Standoff

Therapeutic and Focal Volumes




Front View

*Therapeutic Volume:*

Contour of at least 5 MPa

$$V \approx 2,000 \text{ mm}^3$$

18mm x 20mm x 60 mm




Top View

*Focal Volume:*

Contour of at least 1/2 of Pmax

$$V \approx 10 \text{ mm}^3$$

1.5mm x 2mm x 7 mm



Side View







# CS-Pro VET

10mm Standoff

Therapeutic and Focal Volumes




Front View

*Therapeutic Volume:*

Contour of at least 5 MPa

$$V \approx 1,000 \text{ mm}^3$$

14mm x 20mm x 30 mm




Top View

*Focal Volume:*

Contour of at least 1/2 of Pmax

$$V \approx 10 \text{ mm}^3$$

1.5mm x 2mm x 7 mm



Side View

